

frido

Foldable Travel And Travel⁺ Wheelchair

Designed for effortless mobility



Tool-free assembly

User Manual

DIRECTORY

1. Preface

2. Product Info & Specifications

3. Included Components

4. General Safety Guidelines

5. Assembly & Unboxing

5.1 Unboxing

5.2 Setting up the Wheelchair

6. How to Use the Brake Function

7. How to Fold the Wheelchair

8. Care & Maintenance

8.1 Wheelchair Care

8.2 Maintenance

8.3 Storage & Transportation

9. Warranty

PREFACE

Thank you for choosing the
**Frido Foldable Travel/
Travel⁺ Wheelchair AP-FNN01/AP-FNN02.**

This wheelchair is thoughtfully designed to support seniors, adults, pregnant women, and individuals with various bodily discomforts, disabilities, osteoporosis, post-surgery recovery, and more. It offers a comfortable and convenient mobility solution for those needing assistance with short-distance travel, enhancing daily movement and ease.

- This manual provides important information on setup, safe operation, and maintenance.
- Please read through it carefully to ensure proper use and safety.
- Familiarizing yourself with the features will help you get the most from your wheelchair.
- Our customer support team is always available if you need assistance.

We appreciate your trust in our brand and wish you a safe, enjoyable experience.

About frido

Frido is dedicated to enhancing daily living through innovative, ergonomic solutions that promote comfort and well-being. Our mission is to design the next generation of products to empower people to walk, sit, and sleep better.

Frido's story began with a simple yet powerful mission: to restore dignity and freedom to those who need it most. It all started with a groundbreaking idea—a wheelchair designed not just for mobility but for accessibility and independence. This became the Frido Go, the first-ever self-propelled shower and commode wheelchair capable of navigating complex spaces and folding neatly into a suitcase. From homes to hospitals and even the Indian Army Veterans Rehabilitation Center, this innovation made waves across the world, reaching over eight countries.

But our journey didn't stop there. As we listened to stories of people struggling with daily discomfort, a bigger vision took shape: to create solutions that bring comfort and freedom to everyone, every day. We ventured into designing products that address the smallest pain points and make a big difference—starting from the soles of your feet to the way you sleep at night.

Today, Frido is a trusted name for ergonomic, research-backed products that enhance daily living. Each product we create is fueled by the same passion that drove us to innovate from day one: to give people the freedom to do more. Frido is dedicated to enhancing activities of daily living—be it toileting, bathing, sleeping, or walking—by focusing on human body ergonomics. Our products are thoughtfully designed to alleviate pain and prioritize comfort, redefining how individuals experience everyday tasks.

Frido is a significant player in the consumer healthcare space, Frido is designed in India and made for the world, blending innovation with a global vision to improve lives everywhere.

Product Information & Specifications

PRODUCT	SPECIFICATIONS
Model Name	Frido Foldable Travel and Travel ⁺ Wheelchair
Model Number	AP-FNN01/AP-FNN02
Color	Black
Net Weight	9.5 Kilograms/10.5 Kilograms
Gross Weight	11 Kilograms/13.4 Kilograms
Product Material	Aluminium Alloy
Maximum Climbing Angle	5° slopes
Product Dimensions	32.52L x 19.02W x 35H Inches / 34.5L x 24.5W x 37H Inches (82.6L x 48.3W x 88.9H Centimeters) / (87.6L x 62.2W x 93.9H Centimeters)
Packaging Type	Corrugated Box
Packaging dimensions	24.02L x 15.75W x 31.50H Inches / 25L x 14W x 32H Inches (61.00L x 40W x 80H Centimeters) / (63.5L x 35.5W x 81.2H Centimeters)
Seat Dimensions	~15D X 14W Inches/~15D X 19W Inches
Castor (Front Wheel)	6 inches
Rear Wheel	11 Inches/11.5 Inches

Included Components

ITEM	QUANTITY
Wheelchair	1
User Manual	1
Warranty card	1

4. General Safety Guidelines

1. For the wheelchair to be stable, the rear wheel connecting lever must be placed in the backmost position. Moving the connecting lever can affect how stable and well the wheelchair moves. To stay safe, users should test the wheelchair's stability before using it.
2. Don't use the brakes erratically; it can cause a jerk to the person sitting on the wheelchair.
3. Do not attempt to go up or down stairs while seated on the wheelchair.
4. When sitting in a wheelchair, be careful when going for items outside of the reach of the wheelchair. Avoid leaning too far to prevent overbalancing the wheelchair, which could affect its stability and safety.
5. Before using the wheelchair, please check if the handles or other components are loose. If they are, contact a professional technician for assistance.
6. Do not bend over to pick up the items on the ground.
7. The wheelchair should not be used on the escalator; it can cause serious personal injury.
8. When transferring from the wheelchair to another place, first position the wheelchair properly, with the front wheels facing the direction of transfer. Ensure the wheelchair is parked securely, and check that the brakes are applied before attempting to board.

9. Please do not step on the footrest to get on and off. Be sure to fold the footrest first.

10. When navigating inclines or slopes, avoid abruptly tipping the wheelchair or changing direction quickly. When going downhill, do not brake suddenly, as this can cause the wheelchair to topple forward and create a safety hazard.

11. The maximum weight the wheelchair can bear is 100 kg /120 kg.

12. Wheelchairs are designed for use on flat, even surfaces and should not be driven on uneven terrain. If you encounter an obstacle, navigate around it to avoid tipping the wheelchair.

13. When driving the wheelchair, maintain a constant speed of 3–5 km/h for safe operation.

14. The safe driving range for a wheelchair is a slope between 1.8° and 5°. For steeper slopes, assistance is required. When driving on a slanted surface, the user should be in good physical condition to ensure safe operation.

5. Assembly & Unboxing

5.1 Unboxing

Unpack the folded wheelchair body and user manual from the packaging.

5.2 Setting up the Wheelchair

Step 1: Hold the footrest and the back frame of the wheelchair and push them apart.



Step 2: Open the footrest as shown below.



Step 3: Once open, secure the frame by engaging the side lock/safety lock shown in the image below.



Step 4: Rotate the armrests to their specific position.



6. How to Use the Brake Function



To brake the wheelchair, first pull the brake levers at the back of the chair and then push it upwards until you hear a "snap," indicating the brake is engaged. (shown in the image above)



To release the brake, press the red brake lock cap until you hear another "snap," and the wheelchair will be free to move. (shown in the image above)

7. How to Fold the Wheelchair

1. To fold the wheelchair, rotate the armrests on both sides 270 degrees downward.
2. Then, open the side lock located on the frame.
3. Grip the seat and backrest, and pull them together to fold the wheelchair.
4. Finally, extend the holding handle like a suitcase handle to easily move the chair on its wheels.



8. Care & Maintenance

8.1 Wheelchair Care

1. When cleaning, use a cloth with diluted neutral detergent to wipe and scrub the surface. Then, use a cloth soaked in water to remove any detergent residue and wipe with a dry cloth.
2. Avoid using volatile substances such as paint thinners, gasoline, or other strong chemicals, as they can damage or deteriorate the material.
3. Beware of using disinfectants, as some may corrode metal parts or cause discoloration and deterioration of plastic components.
4. Avoid using ozone or high-pressure sterilizers.

8.2 Maintenance

1. Regularly inspect the wheelchair for safety, ideally every six months, to ensure all connections are secure and the wheelchair operates correctly.
2. Apply appropriate lubricants to frequently moving parts during your regular inspection. Replace parts if necessary.
3. Avoid collisions with the wheelchair or its coating to prevent damage.
4. Keep the wheelchair and its back pads dry, as moisture can cause rust, unusual sounds, or bacterial growth. Wipe away any water immediately.

8.3 Storage & Transportation

8.3.1 Product Storage & Usage Conditions

The manual wheelchair should be stored in a dry and ventilated place and separated from chemically corrosive items such as acid and alkali.

To avoid damage or malfunction, do not place or use the wheelchair in the following conditions:

- Direct sunlight,
- Near heaters,
- Areas with excessive moisture or oil vapors,
- High-temperature, humid, or low-temperature environments,
- Environments with dust, smoke, salt, sulfur, or corrosive substances,
- Unstable or uneven ground.

Ensure the wheelchair is used within a temperature range of 0°C to 40°C and a relative humidity of 80%.

8.3.2 Transportation Conditions

When transporting the wheelchair, always ensure it is properly folded or secured to prevent damage. Use caution when lifting or maneuvering the wheelchair, especially in tight spaces or uneven surfaces.

9. Warranty

From the date of purchase, our company offers repairs and replacements of all wheelchair parts except the wheels to the original purchaser.

Manual Disclaimer:

The images are for illustration purposes only; the actual product may vary. Frido Mobility reserves the sole right to interpret this manual.

Frido



Connect with us

